


Sampling of shrimp pleopods for genotyping on ShrimpID

STEP 1:

Ensure that all equipment, tweezers and scissors, are sterilized prior to use. For information on how to sterilize sampling equipment, view the educational videos on our website www.genics.com or click [here](#).



STEP 2:

Using the scissors, remove a pleopod from the animal.

NOTE: if using broodstock, take pleopods from the back of the animal, closer to the tail.

STEP 3:

Using the scissors, cut the pleopod longitudinally to expose the tissue.

STEP 4:

Place all pleopod pieces in a sample tube with 70% laboratory grade (undenatured) ethanol.

Aim for a 1:1 ratio of tissue to ethanol.

STEP 5:

Record sample vial barcode and sample information (such as identifier tag, tank number, etc.) on the datasheet. Once all sampling has been completed and datasheet has been filled out, prepare samples for shipment to Genics. Shipment can be arranged on the Portal at www.my.genics.com.au.

Contact us

info@genics.com
www.genics.com

Dr. Melony Sellars – CEO

+61 437 025 821
melony.sellars@genics.com.au

Dirk Krueger – CCO

+61 411 019 217
dirk.krueger@genics.com.au

Genics
Be Empowered with your Data